

Laboratory Procedures and Tasks

Table of Contents

Introduction	3
Requisition Slips	3
Patient Information–Requisition Slips	4
Billing–Requisition Slips	5
ABN–Requisition Slips	8
Laboratory–Requisition Slips.....	9
Form Distribution	10
Urine Collection and Phlebotomy	11
Urine Specimen Collection	11
Phlebotomy.....	11
Tubes for Phlebotomy	12
Venipuncture Procedure.....	13
Venipuncture Care	14
Specimen Processing	15
Laboratory and Phlebotomy Area Disinfection.....	17
Laboratory and Bathroom Clean-up	18
Call Quest Diagnostics	18
Maintenance Records	19
Refrigerator Temperature Log	20
Centrifuge Log	21
Bleach Solution	22
References.....	23

Introduction

The purpose of pre-diagnostic cardiovascular screening is to uncover the earliest stages of Cardiovascular Disease (CVD). Current technological advances make it possible to catch this devastating disease when the vessels begin to fail or stiffen. Before these advances, diagnosing **high blood pressure** constitutes the earliest clue that the cardiovascular system suffers damage from inflammation.

Diagnosing CVD at the earliest stages allows early intervention, thereby slowing its advance or even reversing the damage. Simple diet changes like eating more vegetables, limiting sugar and processed foods, promote a healthier body. Including more aerobic activity, reducing stress, and incorporating better eating habits also reduces CVD risk.

Requisition Slips

Specimen collection occurs after a positive Cardio101 test score is above two points. Each patient will have **three** requisition forms, one for each aliquot of specimen (plasma, serum, and urine). See Figure 1.

The image shows a requisition form with several sections and annotations. The form is divided into several sections: 'Bill to section' (top left), 'Fasting box' (middle left), 'Health care provider' (bottom left), 'Patient Information section' (top right), 'Insurance information' (middle right), and 'Test menu' (bottom right). The form contains various fields for patient data, insurance details, and a list of tests.

Annotations on the left side of the form:

- Bill to section
- Fasting box
- Health care provider

Annotations on the right side of the form:

- Patient Information section
- Insurance information
- Test menu

Figure 1. Requisition Form

Patient Information–Requisition Slips

The Medical Laboratory Scientist or Medical Laboratory Technician (MLS/MLT) escorts the patient to the phlebotomy area. Then the MLS/ MLT completes the following tasks:

1. Ask the patient to say their name and birthdate while looking at the patient's license or state ID.
2. Fill out the top, right-hand section of all three-requisition forms with the patient's:

- Last Name, First Name
- Date of Birth
- **Clinic** Phone Number (Figure 2)



Make copies of both the **front** and **back** of the insurance and ID cards.



Always compare the information written on the requisition form with the patient's **license** or **ID**.

Figure 2. Patient Section

3. Verify how many hours the patient fasted. **** Needs to be ≥ 12 hours. ****
4. Mark the appropriate box on the patient's requisition sheet. See Figure 3.
 - If fasting is less than 12 hours – they may need to reschedule the test
 - If the patient brought blood work results done within 6-months, use those results

Fasting Box

Note: Medicare patients may only have certain tests run every six months. An example would be the cholesterol test. Otherwise, the patient pays out of pocket for tests outside of these parameters.

Figure 3. Requisition Form – Fasting Section

Billing–Requisition Slips

Patients will have four different medial insurance scenarios: Private, Medicare, Medicare with a Supplement, and Medicaid.

Private and Senior Select Insurance

1. Fill in the **Billing Section** located at the top, center of the form. See Figure 4.
2. Mark “My Account” for the following:

- Private Insurance
- Senior Select (Care)

Note: Senior Select is a Medicare substitute but their card doesn't say “Medicare.” It is not a supplemental insurance.

Figure 4. Billing Section

Medicare with a Supplement

1. Mark “Medicare” and “Other Insurance” if applicable. See Figure 5.
2. Fill in the Insurance section if not already completed. See Figure 6.

Figure 6. Insurance

Mark the type of insurance

Choose if the patient represents the subscriber, spouse, or other dependent

Patient fills out the ABN form. See page 8.

Figure 5. Billing Section

Billing–Requisition Slips

When patients have Medicare or Medicare with a supplement insurance, they need to fill out an Advanced Beneficiary Notice of Noncoverage (ABN) form. This form states that if test costs are not covered by Medicare or their insurance supplement, the patient will pay the balance due. See page 8.

Medicare

1. Mark “Medicare” in the Billing Section. See Figure 7.
2. Fill in the Insurance section if not already completed. See Figure 8.

Figure 7. Billing Section

Figure 8. Insurance Section

Mark
Insurance
Type

Choose if the
patient represents
the subscriber,
spouse, or other
dependent

Supply
Insurance
Information



Patient fills out an **ABN** form.
See page 8.

Billing–Requisition Slips *Continued*

Patients with Medicaid insurance have requisition forms similar to those with private insurance. **No ABN form is needed.**

Medicaid

3. Mark “Medicaid” in the Billing Section. See Figure 9.
4. Fill out the Insurance section. See Figure 10.

Note: Medicaid insurance usually has “MNSure” after the provider’s name.

Mark
Medicaid

Figure 9. Billing Section

Mark Insurance Type

Choose if the patient
represents the subscriber,
spouse, or other dependent

Supply Insurance Information

Figure 10. Insurance Section

ABN-Requisition Slips

The Advanced Beneficiary Notice of Noncoverage or ABN form is used with the Medicare and Medicare plus supplemental insurance patients. Legally, this document notifies patients that their insurance may not cover certain tests and that they will be responsible for paying the costs. Fill in the laboratory section with the following information for each patient: **Lipid Panel 7600, Glucose 483, and write in (other section) BNP 37386**. See Figure 11.

1. Always filled out with **Medicare Insurance patients**.
2. Photocopy form and keep it with the patient's send out records.
3. Yellow copy- patient keeps for their records.
4. White copy- sent with lab requisitions to Quest.
5. Send a photocopy of patient's ID and insurance cards - front and back to Quest.



If patients have supplemental insurance, the insurance may pick up the cost of the BNP test. Diabetics are only allowed glucose tests quarterly and lipid profiles every 6 months.

<input type="checkbox"/> Hemoglobin A1C <input type="checkbox"/> Lipid Panel <input type="checkbox"/> HbA1c <input type="checkbox"/> HbA1c <input type="checkbox"/> HbA1c	<input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose	<input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose	<input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose <input type="checkbox"/> Fasting Glucose
---	--	--	--

Reason Medicare May Not Pay:

Estimated Cost:

WHAT YOU NEED TO DO NOW:

- Read this notice, so you can make an informed decision about your care.
- Ask us any questions that you may have after you finish reading.
- Choose an option below about whether to receive the checked items listed in the first box above.

Note: If you choose Option 1 or 2, we may help you to use any other insurance that you might have, but Medicare cannot require us to do this.

Options: Check only one box. We cannot choose a box for you.

☒ **OPTION 1:** I want the laboratory test(s) listed above. You may ask to be paid now, but I later want Medicare to pay for an official decision on payment, which is sent to me on a Medicare Summary Notice (MSN). I understand that if Medicare doesn't pay, I am responsible for payment, but I can appeal to Medicare by following the directions on the MSN. If Medicare does pay, you will refund any payments I made to you, less copays or deductibles.

☐ **OPTION 2:** I want the laboratory test(s) listed above, but do not bill Medicare. You may ask to be paid now as I am responsible for payment. I cannot appeal if Medicare is not billed.

☐ **OPTION 3:** I don't want the laboratory test(s) listed above. I understand with this choice I am not responsible for payment, and I cannot appeal to see if Medicare would pay.

Additional Information:

This notice gives our opinion, not an official Medicare decision. If you have other questions on this notice or Medicare billing, call 1-800-MEDICARE (1-800-633-4227/TTY: 1-877-86-2048).

Signing below means that you have received and understand this notice. You also receive a copy.

Signature: Cynthia Snow Date: 6-25-15

Laboratory personnel fill in this section for the Lipid Panel 7600, Glucose 483, and the BNP 37386

Make sure the patient has chosen option 1, 2, or 3

Patients must sign and date the form

Figure 11. ABN Form

Laboratory—Requisition Slips

Because there are three different specimen aliquots, each specimen type is packaged with its own requisition slip. One reason deals with specimen storage. The **Brain Natriuretic Peptide (BNP)** plasma specimen is frozen for transport. The **Lipid Panel**, **Cardio-CRP**, and **Glucose** serum specimens are stored either refrigerated or at room temperature, and the urine **Micro Albumin** test is sent at room temperature.

The largest section on the **Requisition Slip** contains the laboratory test choices. For most patients, the following tests will be marked every time. The exception includes patients who have valid laboratory results from their doctors. See Figures 12 and 13.

Send _____ Client # OR NAME: _____
 Duplicate ADDRESS: _____
 Report to CITY: _____

STATE _____ ZIP _____

PANEL COMPONENTS ON BACK

ORGAN / DISEASE PANELS

34392 ☐ Electrolyte Panel S
 10256 ☐ Hepatic Function Panel S
 10165 ☐ Basic Metabolic Panel w/eGFR S
 10231 ☐ Comp Metabolic Panel w/eGFR S
 7600 ☒ Lipid Panel (Fasting Specimen) S
 14852 ☐ Lipid Panel w/Reflex d-LDL S
 20210 ☐ Obstetric Panel w/Reflex S
 10306 ☐ Hepatitis Panel, Acute w/Reflex S
 10314 ☐ Renal Functional Panel w/eGFR S

HEMATOLOGY

510 ☐ Hemoglobin L
 509 ☐ Hematocrit L
 1759 ☐ CBC (Hgb, Hct, RBC, WBC, Plt) L
 6399 ☐ CBC w/Diff (Hgb, Hct, RBC, WBC, Plt, Diff) L
 8847 ☐ PT with INR L
 763 ☐ PTT, Activated B

OTHER TESTS

7788 ☐ ABO Group & Rh Type Y
 237 ☐ AFP Tumor Marker Y
 779 ☐ _____ Y

4420 ☐ C-Reactive Protein CRP
 29493 ☐ CA 2729
 29256 ☐ CA 125
 303 ☐ Calcium
 10124 ☒ Cardio CRP
 11173 ☐ CCP Ab IgG
 978 ☐ CEA
 334 ☐ Cholesterol, Total
 374 ☐ CK, Total
 375 ☐ Creatinine w/eGFR
 402 ☐ DHEA Sulfate, Immunoassay
 8293 ☐ Estradiol
 4021 ☐ Estradiol
 457 ☐ Ferritin
 466 ☐ Folic Acid
 470 ☐ FSH
 482 ☐ GGT
 8477 ☐ Glucose, Gest. Scr
 8484 ☐ Glucose, Plasma
 483 ☒ Glucose, Serum
 14839 ☐ H. pylori Urea Breath Test
 8435 ☐ hCG, Serum, Qual
 8396 ☐ hCG, Serum, Qual

Fill in the boxes

10124 Cardio CRP
(also known
as hs-CRP)

600 Lipid Panel

483 Glucose

ZIP _____

C-Reactive Protein CRP S
 CA 2729 S
 CA 125 S
 Calcium S
 Cardio CRP S
 CCP Ab IgG S
 CEA S
 Cholesterol, Total S
 CK, Total S
 Creatinine w/eGFR S
 DHEA Sulfate, Immunoassay S
 Direct LDL S
 Estradiol S
 Ferritin S
 Folic Acid S
 FSH S
 GGT S
 Glucose, Gest. Scr GY
 Glucose, Plasma GY
 Glucose, Serum S
 H. pylori Urea Breath Test HB
 hCG, Serum, Qual S
 hCG, Serum, Qual L
 Hemoglobin A1c L

571 ☐ Iron, Total S
 593 ☐ LDH S
 599 ☐ Lead (B) TN
 615 ☐ Lf S
 606 ☐ Lipase S
 6646 ☐ Lyme Disease Ab w/Reflex to Bld IgG S
 622 ☐ Magnesium S
 6517 ☒ Microalbumin, Random Urine w/Creat
 11290 ☐ DX
 11293 ☐ MCR Scr
 718 ☐ Phosphorus S
 733 ☐ Potassium S
 745 ☐ Progesterone S
 746 ☐ Prolactin S
 5363 ☐ PSA, Total S
 793 ☐ Reticulocyte Count, Automated L
 4418 ☐ Rheumatoid Factor S
 799 ☐ RPR (Monitoring) w/Reflex Titer S
 36125 ☐ RPR (DX) w/Reflex Confirm S
 802 ☐ Rubella IgG S
 809 ☐ Sed Rate by Mod West L
 15983 ☐ Testosterone, Total, LC/MS/MS SR
 873 ☐ Testosterone, Total, Male SR

Fill in the box

6517 Microalbumin,
Random Urine
with Creat.

Figure 12. Requisition Form Number One

Figure 13. Requisition Form Number Two

Laboratory–Requisition Slips *Continued*

The third requisition slip will be the same for every patient. Write the following in the space located beneath the lab test section. See Figure 14.

3 7386 - BNP - frozen Specimen

37386 – BNP – Frozen Specimen

Figure 14. Requisition Slip Number Three

Form Distribution

At this point, there are three requisition slips, one ABN form if applicable, and two copies of the patient's insurance and ID cards. After the patient fills out the ABN and the requisition slips are completed, photocopy them too. The following table outlines where each of these documents go after the specimens are obtained. See Table 1.

Table 1. Document Distribution

Documents	Patient	Patient's Chart	Lab Record Book	Specimen
Requisition Slip				One for each specimen type (3)
Requisition Slip photocopy			✓	
ABN form	Yellow copy only			White copy
ABN form photocopy			✓	
Insurance and Photo ID copy		✓		✓ (Medicare insurance)

Urine Collection and Phlebotomy

Urine Specimen Collection

The microalbumin creatinine ratio tests kidney function. Albumin is a protein found in the blood. Normally, no albumin should be excreted by patients with normal kidney function. This test is usually done on patients who have diabetes, high blood pressure, or kidney disease. Creatinine excretion happens at a steady rate throughout the day, while albumin amounts vary.

1. Hand the patient a urine specimen cup labeled with their name, birthdate, date, and time collected. See Figure 15.
2. Escort them to the bathroom.
3. Show the patient where the specimen pass-through door is located.
4. Invite the patient to sit in the phlebotomy chair when finished.



Figure 15. Urine Specimen Cup
<https://www.youtube.com/watch?v=T9eOPhpUeF8&t=53s>

Phlebotomy

Phlebotomy is the procedure in which blood samples are drawn through venipuncture. Four categories of laboratory professionals are proficient in venipuncture. They are Medical Assistants (MAs), Phlebotomists, Medical Laboratory Technicians (MLTs), and Medical Laboratory Scientists/ Medical Technologists (MLS/MT). Sometimes the Physician Assistant (PA) draws blood samples.

It is important to gather all supplies needed for the procedure. Check the phlebotomy tray for needed supplies and replace them before the first scheduled patient arrives. Always have extra tubes available in case there is an issue. Examples of issues are an incomplete blood draw or the patient's blood clots before the tube is filled.



Phlebotomy *Continued*

Tubes for Phlebotomy

The only three tubes will be drawn on each patient. **Two SST Tubes** or red and gray mottled gel tubes are used for the glucose, lipid panel, and the hs-CRP tests. The SST or serum separator tubes will clot for about 30 to 60 minutes, spin in the centrifuge for 15 minutes, and then is bagged for transport. They do not need to be poured over into a transport tube.

One EDTA Lavender-Top Tube is drawn for the BNP test. Only one lavender-top tube is needed to collect enough plasma. The lavender-top tube is spun as soon as possible, the plasma pipetted off, and put into a screw-top transport container. Place the labeled specimen in the freezer. See Table 2. For specimen processing, see page 15.

Table 2. Laboratory Draw Needs

Number of Tubes	Tube Type	Tests	Sample Needed	Storage
	SST Tube	Lipid Panel	Serum	Room Temperature
	SST Tube	Glucose, hs-CRP	Serum	Room Temperature
	Lavender – Top Tube	BNP	Plasma	Freeze

Phlebotomy Continued

Venipuncture Procedure

1. Ask the patient to state their name and birthdate.
2. Verify information matches the patient requisitions and labels.
3. Prepare all supplies needed for venipuncture.
4. Ask the patient to roll up their sleeve to access the veins.
5. Talk about the patient's pet or favorite vacation spot to help them relax.
6. Put gloves on.

Blood Specimen Collection Process

1. Tie the tourniquet around the upper arm.
2. Wipe antecubital area with an alcohol wipe using concentric circles/spirals.
3. Uncap needle while the alcohol dries.
4. Draw **2 SST tubes** (red-topped tubes with gel in the bottom).
5. Draw **1 full lavender-top tube**.
6. Mix tubes by gentle inversion 10 times.

Tubes Filled

1. Release the tourniquet.
2. Grab gauze squares with a free hand and place gently over puncture site.
3. Remove needle quickly, then apply pressure over the venipuncture site.

Labeling tubes

1. Ask the patient to hold pressure on the puncture site.
2. Write phlebotomist's initials, date, and time drawn on the labels.
3. Place labels on the tubes.
4. Ask the patient to read out loud the name and birthdate on the tubes.
5. Place specimens in the phlebotomy tray

Venipuncture Care

1. Lift up the gauze to access bleeding, then replace the gauze and continue holding pressure on the site.
2. Take a piece of paper tape and place it over the site if bleeding has stopped.
3. Instruct the patient to keep pressure on the site for 5 minutes.
4. Escort the patient to the Vascular Technician (VT) or to the 4-Test Procedure.



If bleeding doesn't stop, apply pressure for a few more minutes. Recheck for bleeding. Repeat steps 1 through 3 until bleeding stops.

When patients are on anticoagulant therapy, their blood may take longer to stop bleeding. The other common reasons for continued patient bleeding are:

- Arthritis medication or medicine containing aspirin
- Low platelet count
- Accidentally drawing an artery



Remember that arteries have a pulse. **Do not draw an artery.**

Specimen Processing

Specimen integrity may be affected when plasma or serum remains in contact too long with RBC's or red blood cells. The BNP or brain natriuretic peptide results tend to rise with heart failure. Therefore, it is important to follow the processing guidelines to ensure accurate test results. See Table 3.

1. Centrifuge **lavender-top tubes** within 15 minutes of the blood draw.
2. Separate plasma from the RBCs within 30 minutes of the blood draw or test results will be affected.
3. Clotting in **SST tubes** takes 30 to 60 minutes.
4. Tilt **SST tubes** to see if the serum in has separated from the clot.
5. Centrifuge **SST tubes** when ready.

Table 3. Specimen Processing

Test Name	When to Centrifuge	Minutes for Centrifuging	Specimen type: Plasma, Serum or Urine	Storage Requirements	Medical Condition
Lipid Profile Total Cholesterol HDL LDL Total Chole/HDL ratio Triglycerides % of chylomicrons, VLDL, LDL, HDL	30 to 60 minutes and the specimen are clotted	15 minutes	Serum	Room Temperature	CVD – cardiovascular disease
Glucose	30 to 60 minutes until clotted	15 minutes	Serum	Room Temperature	CVD/ Diabetes
Hs-CRP	30 to 60 minutes until clotted	15 minutes	Serum	Room Temperature	Shows persistent levels of inflammation—CVD
BNP	Spin within 15 minutes	15 minutes	Plasma	Frozen	CVD
Micro-albumin with Creatinine Ratio	N/A	N/A	Urine	Room Temperature	Kidney disease

Specimen Processing *Continued*

Lavender Top Tube

1. Mark transfer vial *plasma*, and place specimen and requisition labels on the vial.
2. Put on PPE (personal protective equipment: disposable lab coat, gloves, and face shield).
3. Pipette the plasma off the RBCs in the **lavender-topped tube** and transfer to the labeled vial.
4. Place capped specimen in the large pocket of the transport bag marked *frozen specimen* (orange stickers).
5. Make sure specimen label barcode is visible through the pouch.
6. Photocopy all requisition and ABN slips, mark **Lab Test Sent** on copies and place in 3-ring binder marked *Lab Test Sent*.
 - Place original corresponding paperwork in the pocket (barcode is visible).
 - Lab **ABN** form (if needed, see appendix)
 - **Requisition Slip**
7. Place specimen in the freezer ASAP (as soon as possible).

SST Tubes

1. Check **SST tubes** for clotting (serum should be visible when the tube is tilted).
2. Place specimen label from the requisition on each tube before spinning.
3. Spin tubes for 15 minutes in the centrifuge.
4. Remove 2 SST tubes from one patient at a time.
5. Place both tubes into one transport bag labeled *room temperature* (large pocket).
 - Photocopy all requisition and ABN slips, mark **Lab Test Sent** on copies and place in 3-ring binder marked *Lab Test Sent*
 - Place original corresponding paperwork in the pocket (barcode is visible).
 - Lab **ABN** form (if needed: see appendix)
 - **Requisition Slip**
 - Leave at room temperature until specimen pick-up

Specimen Processing *Continued*

Urine

1. Mix **urine** specimen by gently swirling collection cup about 3 to 5 times.
2. Label transport urine vial with the patient ID and requisition label.
3. Pipette urine and fill transport vial to the line and screw on the cap.
4. Place vial in transport bag labeled *room temperature* (large pocket).
 - Photocopy all requisition and ABN slips, mark **Lab Test Sent** on copies and place in 3-ring binder marked *Lab Test Sent*
 - Place original corresponding paperwork in the pocket (barcode is visible)
 - Lab **ABN** form (if needed: see appendix)
 - **Requisition Slip**
 - Leave at room temperature until specimen pick-up

Laboratory and Phlebotomy Area Disinfection

While tubes are spinning and clotting, cleaning and disinfection can be done in the phlebotomy area.

1. Disinfect phlebotomy area (if last draw is completed or if blood and body fluids spilled).
2. Carry needle waste container back to the lab.
3. Disinfect the phlebotomy chair with CaviWipes XL.
4. Disinfect countertops and floor area with a 1:16 bleach solution (made fresh each morning).



See page 21 for making the 1:16 bleach solution.

Laboratory and Bathroom Clean-up

1. Empty urine specimens in the toilet (at shift end).
2. Destroy labels and discard urine cups in the appropriate waste container.
3. Disinfect bathroom:
 - Use CaviWipes XL on specimen pass through, light switch, faucet handles, and door knob
 - Use 1:16 bleach solution for the toilet, the floor around the toilet, and sink
4. Disinfect Lab.
 - Use CaviWipes XL on stainless steel work benches and centrifuge (inspect rotors for cracks)
 - Use 1:16 bleach solution on the floor
5. Disinfect specimen pick-up areas on the bench and refrigerator with CaviWipes XL after specimen pick-up.
6. Turn off lights and shut the lab door.

Call Quest Diagnostics

1. Call Quest Diagnostics for specimen pick-up **before 11:00 am**. 1-866-697-8378.
 - Select option1 (for healthcare providers)
 - Select option 6 for specimen pick-up (have clinic ID number ready)



See instructions for making fresh 1:16 bleach solution daily on page 21.

Maintenance Records

Refrigerator

Recording the refrigerator and freezer temperatures daily ensures proper temperatures for specimen storage. **Frozen** specimens, such as the plasma, need to be maintained between **0° to 14° F**, and the **refrigerator** temperature is between **35° to 45° F**. Record the temperatures on the log sheet. See page 19.

Centrifuge

When properly maintained, used, and inspected, the **centrifuge** remains a safe piece of equipment. However, severe accidents may occur if the centrifuge is allowed to run while unbalanced. Running a centrifuge while unbalanced causes the rotors to wear and break, damaging the centrifuge. It is important to consult and follow the manufacturers instructions for routine checks and care.

After **500** hours of use, the Octofuge will be sent to the manufacturer and be inspected by a trained technician. They will inspect and test the centrifuge for any defects, worn bearings, and rotors. The speed will be checked with a tachometer to ensure the speed matches 3,000 rpms. See Figure 16.

See page 20 for the log sheet and record the date, hours, and maintenance done.




Figure 16. Example of a Six Place Swing Out Rotor Without Buckets
Source: <https://www.youtube.com/watch?v=LrBsS1agd2A>

Refrigerator Temperature Log

[illegible]

Centrifuge Log

[illegible]

Bleach Solution

The following **1:16 bleach solution** is considered an intermediate-strength germicide. According to Scungio (2017), this bleach solution "will eliminate most bacteria (including *Mycobacterium tuberculosis*) and all fungi, and it inactivates viruses" (Routine Decontamination, para 2). To ensure proper efficacy, the bleach solution is made **daily**.

1. Use protective equipment when diluting bleach.
 - Lab coat
 - Gloves
 - Chemistry goggles
 - Hood with an exhaust fan
2. Gather supplies
 - Bleach with **8.25%** sodium hypochlorite
 - Distilled water
 - Graduated cylinders (1-500 ml and 1-100 ml)
 - Funnel
 - Laboratory safety wash bottle (16 oz or 500 ml)
3. Measure **444 ml of distilled water** using the **500 ml graduated cylinder**.
4. **Pour 244 ml of distilled water** into the laboratory **safety wash bottle** using the funnel.
5. Measure **30 ml of 8.25% bleach** into the **100 ml graduated cylinder** and pour into the safety wash bottle.
6. Pour in the remaining **200 ml of distilled water** into the laboratory safety wash bottle, cap and gently swirl to mix.
7. Label the bottle with the date made.

443 ml distilled water + 30 ml of 8.25% bleach = 1:16 dilution of 8.25% bleach solution



Remember to add **acids** to **water**.

References

- Cohn, Jay N., Duprez, Daniel A., Grandits, Gregory A. (2005, July). Arterial elasticity as part of a comprehensive assessment of cardiovascular risk and drug treatment. *Hypertension*, 2017-2020.
doi:10.1161/01.HYP.0000165686.50890.c3
- Miller, Michael J., Astles, Rex, Baszler, Timothy, Chapin, Kimberle, Carey, Roberta . . . Wiedbrauk, Danny. (2012, January 6). Guidelines for safe work practices in human and animal medical diagnostic laboratories: Recommendations of a CDC-convened, biosafety blue ribbon panel. In *Morbidity and Mortality Weekly Report (MMWR)*. Retrieved from https://www.cdc.gov/mmwr/preview/mmwrhtml/su6101a1.htm?s_cid=su6101a1_w
- Scungio, Daniel. (2017, June). Lab decontamination tools and practices. *Medical Laboratory Observer*. Retrieved from <https://www.mlo-online.com/diagnostics/pathology/article/13009199/lab-decontamination-tools-and-practices>